

product number

60/66 Bb

## A-Trisaccharide-Biotin-BSA

GalNAcα1-3Galβ1-O-APE-Biotin-BSA 2 Fucα1

The oligosaccharide is linked to biotinylated-BSA lysine reidues via a spacer using the isothiocyanate coupling method. The spacer used is *p*-aminophenylethyl (APE). The spacer is glycosidically linked to the reducing monosaccharide unit of the oligosaccharide.

The conjugate normally contain 10-20 moles oligosaccharide per mole BSA.

The conjugate is supplied as freeze dried powder. Storage: -18° C.